

## Pharmacy Refrigerators HYC-290



### Ergonomic Design



High Accuracy  
Temperature Control



Interior  
LED



Lockable  
Design



### Air Cooling, Auto Defrost

Better cabinet temperature uniformity,  
auto defrost



### Adjustable Shelves

The shelves are adjustable to any height to  
fulfill different requirements



### USB Interface

Data of temperature history can be downloaded  
via the built-in USB port and optional printer design



### Heated Double Pane Glass Door

Clear display of stored products without condensation  
at 32 °C and 85% relative humidity condition



### Self Closing Door with 90° Stay Open Feature

Energy conserving, maintains temperature  
uniformity to protect valuable samples

# Specifications

## HYC-290

### Technical Data

|               |                            |
|---------------|----------------------------|
| Cabinet Type  | Upright, single glass door |
| Climate Class | N                          |
| Cooling Type  | Forced air cooling         |
| Defrost Mode  | Auto                       |
| Refrigerant   | CFC-Free                   |
| Noise((dB(A)) | 43                         |

### Performance

|               |     |
|---------------|-----|
| Temp Range(℃) | 2~8 |
|---------------|-----|

### Control

|            |                |
|------------|----------------|
| Controller | Microprocessor |
| Display    | LED            |

### Electrical Data

|                       |               |        |
|-----------------------|---------------|--------|
| Power Supply(V/Hz)    | 220-240/50/60 | 115/60 |
| Power(W)              | 340           | 460    |
| Electrical Current(A) | 2.2           | 5.0    |

### Dimensions

|                              |                                     |
|------------------------------|-------------------------------------|
| Capacity(L/Cu.Ft)            | 290/10.2                            |
| Net/Gross Weight(approx)     | 105/127(kg) 231.5/280.0(lbs)        |
| Interior Dimensions(W*D*H)   | 530*555*1080(mm) 20.9*21.9*42.5(in) |
| Exterior Dimensions(W*D*H)   | 665*710*1665(mm) 26.2*28.0*65.6(in) |
| Packing Dimensions(W*D*H)    | 735*800*1800(mm) 28.9*31.5*70.9(in) |
| Container Load(20'/40'/40'H) | 21/42/42                            |



### Functions

|               |   |
|---------------|---|
| High/Low Temp | Y |
| Remote System | Y |
| Power Failure | Y |
| Sensor Error  | Y |
| Low Battery   | Y |
| Door Ajar     | Y |

### Accessories

|                 |     |
|-----------------|-----|
| Caster          | Y   |
| Foot            | Y   |
| Test Hole       | Y   |
| Shelves /Basket | 5/- |
| USB Interface   | Y   |

### Others

|             |       |
|-------------|-------|
| Certificate | CE,UL |
|-------------|-------|

## Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1  
Haier Road, Qingdao, 266101, P.R. China  
Tel: +86-0532-88936011 / 5955  
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

